

UC-331 Bathroom Scale

Specifications;

- Capacity: 150kg
- Measurement Range: 5kg ~ 150kg
- Display Resolution: 100g
- Display Type: Liquid crystal
- Battery: Coin type lithium, CR2032 x 1
- Battery Life: Approximately 2000 readings (with new batteries)
- Operating Temperature: 10 to 40 degrees Celsius
- Storing Temperature: 0 to 50 degrees Celsius
- Dimensions: 310mm x 300mm x 20mm
- Weight: 1.5kg

Before Use;

Please read these instructions before use. The battery is protected from the point of manufacture by insulation tape. Before use, please remove the tape which is exposed from the battery lid on the rear of the scale. The battery supplied with the scale is for initial calibration and set up. A&D cannot guarantee the life of the battery supplied with the scale, and it is not covered by warranty.

Instructions;

1. Place the scale on a solid flat surface; ensure the scale is dry to avoid the potential for slipping. (Do not measure on carpet).
2. Tap the top of the scale to activate weighing mode. The unit will display a series of "8" then "0.0" will appear on the display.
3. Step on the scale, and remain as still as possible.
4. The display will stabilise, and the measurement will be displayed.
5. Step off the scale and the measurement will be locked on for 6 seconds.
6. The scale will automatically turn off.

Error Messages;

"OL" Overload

"LO" Low Battery

"----" No measurement is possible as you stepped on the scale before "0.0" was displayed.

Maintenance;

Clean with a damp cloth, do not pour any liquid directly on the scale. Remove the battery during long periods without use. Store the scale in a clean dry environment.

Warranty;

This scale has a warranty period of 12 months from date of purchase. (Please refer to A&D Australasia's warranty document for details).